


Anti-Grid1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ulk0

Item Number: CM0026807



Introduction

High-quality rabbit polyclonal antibody targeting human GRID1 (GluD1), validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Suitable for studies on synaptic organization and plasticity.

[Learn More](#)

Parameter	Specification
Product Name	GRID1 Rabbit pAb
Remarks/Alias	GluD1; GluD1-b; GRID1
Species	Human
GeneID (Human)	2894
GeneID	2894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 880-1009 of human GRID1 (NP_060021.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9ULK0
Protein Weight	112kDa
Shipping	Ice bag