


Anti-Grid2 Polyclonal Antibody For Wb, If/Icc, Elisa - O43424

Item Number: CM0021872



Introduction

Rabbit polyclonal antibody against human GRID2 (GluD2), suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying synaptic plasticity and cerebellar function.

[Learn More](#)

Parameter	Value
Product Name	GRID2 Rabbit pAb
Remarks/Alias	GluD2; SCAR18; GRID2
Species	Human
GeneID (Human)	2895
GeneID	2895
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 900-1007 of human GRID2 (NP_001501.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43424
Protein Weight	113kDa
Shipping	Ice bag