

Anti-Grik4 Polyclonal Antibody For Wb, Elisa - Q16099

Item Number: CM0018236

Introduction

GRIK4 Rabbit Polyclonal Antibody validated for WB and ELISA; reacts with human, mouse, and rat. Targets the kainate receptor subunit GRIK4 (KA1), suitable for neuroscience research on glutamate signaling.

[Learn More](#)

Parameter	Value
Product Name	GRIK4 Rabbit pAb
Remarks/Alias	KA1; EAA1; GRIK; GluK4; GluK4-2; GRIK4
Species	Human
GeneID (human)	2900
GeneID	2900
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-220 of human GRIK4 (NP_055434.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16099
Protein Weight	107kDa
Shipping	Ice bag