

Anti-Grin2B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13224

Item Number: CM0030460



Introduction

High-quality rabbit polyclonal antibody against GRIN2B (NMDAR2B). Suitable for WB, IF/ICC, ELISA. Reacts with human, mouse, and rat. Ideal for neuroscience research and ion channel studies.

[Learn More](#)

Parameter	Value
Product Name	GRIN2B Rabbit pAb
Remarks/Alias	NR3; MRD6; NR2B; hNR3; DEE27; EIEE27; GluN2B; NMDAR2B; GRIN2B
Species	Human
Gene ID (Human)	2904
Gene ID	2904
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1185-1484 of human GRIN2B (NP_000825.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13224
Protein Weight	166kDa
Shipping	Ice bag