


Anti-Grin2D Rabbit Polyclonal Antibody For Wb, Elisa - O15399

Item Number: CM0018245



Introduction

GRIN2D rabbit polyclonal antibody validated for WB and ELISA, reactive with human and mouse. Ideal for neuroscience research on NMDA receptor function and synaptic plasticity.

[Learn More](#)

Parameter	Specification
Product Name	GRIN2D Rabbit pAb
Remarks/Alias	EB11; NR2D; DEE46; EIEE46; GluN2D; NMDAR2D; GRIN2D
Species	Human
GeneID (human)	2906
GeneID	2906
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 330-480 of human GRIN2D (NP_000827.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15399
Protein Weight	144kDa
Shipping	Ice bag