

Anti-Grina Polyclonal Antibody For Wb, Elisa - Q7Z429

Item Number: CM0022993

Introduction

Rabbit polyclonal antibody against human GRINA (Protein lifeguard 1), validated for WB and ELISA. Detects human, mouse, and rat GRINA. Ideal for apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	GRINA Rabbit pAb
Remarks/Alias	LFG1; HNRGW; TMBIM3; NMDARA1; GRINA
Species	Human
GeneID (Human)	2907
GeneID	2907
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GRINA (NP_000828.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z429
Protein Weight	41kDa
Shipping	Ice bag