

Anti-Grip1 Polyclonal Antibody For Wb, Elisa - Q9Y3R0

Item Number: CM0029309



Introduction

High-quality rabbit anti-GRIP1 polyclonal antibody validated for WB and ELISA, with cross-reactivity in human and mouse. Ideal for studying glutamate receptor scaffolding and neuronal signaling complex assembly.

[Learn More](#)

Parameter	Value
Product Name	GRIP1 Rabbit pAb
Remarks/Alias	GRIP; FRASRS3; GRIP1
Species	Human
GeneID (human)	23426
GeneID	23426
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human GRIP1 (NP_001171545.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y3R0
Protein Weight	122kDa
Shipping	Ice bag