

Anti-Grk1 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q15835

Item Number: CM0026106



Introduction

High-quality rabbit polyclonal antibody targeting human GRK1 (rhodopsin kinase). Validated for WB, IF-P, IHC-P, and ELISA with reactivity across human, mouse, and rat. Ideal for phototransduction and retinal research.

[Learn More](#)

Parameter	Value
Product Name	GRK1 Rabbit pAb
Remarks/Alias	RK; RHOK; GPRK1; GRK1
Species	Human
GeneID (Human)	6011
GeneID	6011
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human GRK1 (NP_002920.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15835
Protein Weight	64kDa
Shipping	Ice bag