

Anti-Grk1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15835

Item Number: CM0029103



Introduction

Anti-GRK1 rabbit polyclonal antibody for Western blot, immunohistochemistry, and ELISA applications. Detects human, mouse, and rat GRK1; immunogen targeting N-terminal region. For research in vision and retinal biology.

[Learn More](#)

Parameter	Value
Product Name	GRK1 Rabbit pAb
Remarks/Alias	RK; RHOK; GPRK1; GRK1
Species	Human
Gene ID (Human)	6011
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human GRK1 (NP_002920.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15835
Protein Weight	64kDa
Shipping	Ice bag