

Anti-Grk2 Monoclonal Antibody [Ko Validated] For Wb, If/Icc, Elisa - P25098

Item Number: CM0011968



Introduction

Rabbit monoclonal anti-GRK2 antibody validated for knockout, suitable for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for GPCR desensitization and beta-arrestin signaling studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GRK2 Rabbit mAb
Remarks/Alias	BARK1, ADRBK1, BETA-ARK1, K2
Species	Human
GeneID	156
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 590-689 of human GRK2 (P25098).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25098
Protein Weight	80 kDa
Shipping	Ice bag