

Anti-Grk2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P25098

Item Number: CM0022124



Introduction

Rabbit polyclonal antibody against human GRK2 (beta-adrenergic receptor kinase 1). Validated for WB, IHC-P, ELISA. Reacts with Human, Mouse, Rat. Ideal for GPCR signaling and desensitization studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GRK2 Rabbit pAb
Remarks/Alias	BARK1; ADRBK1; BETA-ARK1; K2
Species	Human
Gene ID (Human)	156
Gene ID	156
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 440-689 of human GRK2 (NP_001610.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25098
Protein Weight	80kDa
Shipping	Ice bag