

Anti-Grk3 Rabbit Monoclonal Antibody For Wb, Elisa - P35626

Item Number: CM0010556

Introduction

High-quality rabbit monoclonal antibody targeting GRK3 (BARK2). Validated for WB and ELISA with human cross-reactivity. Essential tool for GPCR desensitization and signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	GRK3 Rabbit mAb
Remarks/Alias	BARK2; ADRBK2; GRK3
Species	Human
GenelD (Human)	157
GenelD	157
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P35626
Protein Weight	80kDa
Shipping	Ice bag