

Anti-Grk5 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P34947

Item Number: CM0027242



Introduction

High-quality GRK5 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with Human, Mouse, Rat. Essential tool for studying GPCR kinase 5 function.

[Learn More](#)

Parameter	Value
Product Name	GRK5 Rabbit pAb
Remarks/Alias	GPRK5; FP2025; GRK5
Species	Human
GeneID (Human)	2869
GeneID	2869
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34947
Protein Weight	68kDa
Shipping	Ice bag