

Anti-Grk6 Polyclonal Antibody For Wb, Elisa - P43250

Item Number: CM0028672



Introduction

Rabbit polyclonal antibody targeting human GRK6 (66 kDa). Suitable for WB and ELISA. Cross-reacts with mouse and rat. Used in GPCR desensitization studies. UniProt P43250.

[Learn More](#)

Parameter	Value
Product Name	GRK6 Rabbit pAb
Remarks/Alias	GPRK6; GRK6
Species	Human
GeneID (Human)	2870
GeneID	2870
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 268-465 of human GRK6 (NP_001004106.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43250
Protein Weight	66kDa
Shipping	Ice bag