

Anti-Grm1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13255

Item Number: CM0017533



Introduction

Rabbit polyclonal antibody targeting human GRM1 (metabotropic glutamate receptor 1), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience research on synaptic plasticity and CNS signaling.

[Learn More](#)

Parameter	Value
Product Name	GRM1 Rabbit pAb
Remarks/Alias	MGLU1; SCA44; GPRC1A; MGLUR1; SCAR13; PPP1R85; GRM1
Species	Human
Gene ID (Human)	2911
Gene ID	2911
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 920-1140 of human GRM1 (NP_001264993.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13255
Protein Weight	132kDa
Shipping	Ice bag