

Anti-Grm4 Rabbit Polyclonal Antibody For Wb And Elisa - Q14833

Item Number: CM0028801



Introduction

High-quality rabbit polyclonal antibody targeting human GRM4 (mGlu4, GPRC1D, MGLUR4). Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for metabotropic glutamate receptor 4 research.

[Learn More](#)

Parameter	Value
Product Name	GRM4 Rabbit pAb
Remarks/Alias	mGlu4, GPRC1D, MGLUR4, GRM4
Species	Human
Gene ID	2914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-290 of human GRM4 (NP_000832.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q14833
Protein Weight	102kDa
Shipping	Ice bag