

Anti-Grm5 Rabbit Monoclonal Antibody For Wb, Elisa - P41594

Item Number: CM0005792



Introduction

Highly specific anti-GRM5 rabbit monoclonal antibody for WB and ELISA. Recognizes human, mouse, and rat metabotropic glutamate receptor 5, a key GPCR in synaptic plasticity.

[Learn More](#)

Parameter	Value
Product Name	GRM5 Rabbit mAb
Remarks/Alias	mGlu5, GPRC1E, MGLUR5, PPP1R86, GRM5
Species	Human
Gene ID (Human)	2915
Gene ID	2915
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1100-1200 of human GRM5 (P41594)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P41594
Protein Weight	132kDa
Shipping	Ice bag