

Anti-Grm5 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P41594

Item Number: CM0024741



Introduction

Rabbit polyclonal antibody against human metabotropic glutamate receptor 5 (GRM5/mGlu5). Validated for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience and synaptic plasticity studies.

[Learn More](#)

Parameter	Details
Product Name	GRM5 Rabbit pAb
Remarks/Alias	mGlu5; GPRC1E; MGLUR5; PPP1R86; GRM5
Species	Human
GeneID	2915
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 870-1070 of human GRM5 (NP_001137303.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41594
Protein Weight	132kDa
Shipping	Ice bag