


Anti-Gro Alpha Polyclonal Antibody For Wb, If/Icc, Elisa - P09341

Item Number: CM0030458



Introduction

Rabbit polyclonal antibody targeting human GRO alpha (CXCL1), validated for WB, IF/ICC, and ELISA. Cross-reacts with rat. For inflammation and neutrophil chemotaxis research.

[Learn More](#)

Parameter	Value
Product Name	GRO alpha Rabbit pAb
Remarks/Alias	FSP; GRO1; GROa; MGSA; NAP-3; SCYB1; MGSA-a; GRO alpha
Species	Human
GeneID (Human)	2919
GeneID	2919
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-107 of human GRO alpha (NP_001502.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P09341
Protein Weight	11kDa
Shipping	Ice bag