


Anti-Grp94 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P14625

Item Number: CM0014292



Introduction

Rabbit polyclonal antibody against human Grp94 (Endoplasmic), validated for WB, IF/ICC, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Targets the 92 kDa ER chaperone involved in protein folding and TLR signaling.

[Learn More](#)

Parameter	Value
Product Name	Grp94 Rabbit pAb
Remarks/Alias	ECGP; GP96; TRA1; GRP94; HEL35; HEL-S-125m; Grp94
Species	Human
Gene ID (Human)	7184
Gene ID	7184
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 40-310 of human Grp94 (NP_003290.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14625
Protein Weight	92kDa
Shipping	Ice bag