


Anti-Grp94 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P14625

Item Number: CM0025199



Introduction

Grp94 rabbit polyclonal antibody (SWISS P14625) validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat endoplasmic reticulum chaperone involved in protein folding and tumor immunology.

[Learn More](#)

Parameter	Value
Product Name	Grp94 Rabbit pAb
Remarks/Alias	ECGP, GP96, TRA1, GRP94, HEL35, HEL-S-125m, Grp94
Species	Human
Gene ID (Human)	7184
Gene ID	7184
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 501-750 of human Grp94 (NP_003290.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14625
Protein Weight	92kDa
Shipping	Ice bag