


Anti-Grpel1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Hav7

Item Number: CM0026261



Introduction

Rabbit polyclonal antibody targeting human GRPEL1 (mitochondrial GrpE homolog), suitable for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial protein import and chaperone research.

[Learn More](#)

Parameter	Value
Product Name	GRPEL1 Rabbit pAb
Remarks/Alias	GrpE; HMGE; mt-GrpE#1; GRPEL1
Species	Human
GeneID (Human)	80273
GeneID	80273
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-217 of human GRPEL1 (NP_079472.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HAV7
Protein Weight	24kDa
Shipping	Ice bag