


Anti-Grpr Rabbit Polyclonal Antibody For Wb And Elisa - P30550

Item Number: CM0025318



Introduction

High-quality rabbit polyclonal antibody against human GRPR, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying GPCR signaling, itch sensation, and respiratory control.

[Learn More](#)

Parameter	Value
Product Name	GRPR Rabbit pAb
Remarks/Alias	BB2; BB2R; BRS2; GRPR
Species	Human
Gene ID (Human)	2925
Gene ID	2925
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 235-384 of human GRPR (NP_005305.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30550
Protein Weight	43kDa
Shipping	Ice bag