


Anti-Grsf1 Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q12849

Item Number: CM0002883



Introduction

GRSF1 Rabbit monoclonal antibody validated for WB, IF/ICC, IP, ELISA. Reacts with human, mouse, rat GRSF1, a mitochondrial post-transcriptional gene expression regulator. Ideal for mitochondrial ribosome assembly and RNA binding studies.

[Learn More](#)

Parameter	Value
Product Name	GRSF1 Rabbit mAb
Remarks/Alias	GRSF1
Species	Human
Gene ID (Human)	2926
Gene ID	2926
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human GRSF1 (Q12849)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12849
Protein Weight	53 kDa
Shipping	Ice bag