


Anti-Grsf1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q12849

Item Number: CM0022777



Introduction

Rabbit polyclonal antibody targeting G-rich sequence factor 1 (GRSF1), validated in WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat; essential for mitochondrial RNA biology studies.

[Learn More](#)

Parameter	Value
Product Name	GRSF1 Rabbit pAb
Remarks/Alias	GRSF1
Species	Human
GeneID (human)	2926
GeneID	2926
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human GRSF1 (NP_002083.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12849
Protein Weight	53kDa
Shipping	Ice bag