


Anti-Gsc Monoclonal Antibody For Wb, Elisa - P56915

Item Number: CM0002014



Introduction

Anti-GSC rabbit monoclonal antibody validated for Western blot and ELISA. Targets homeobox protein goosecoid (P56915) in human and mouse, ideal for developmental biology and organogenesis research.

[Learn More](#)

Parameter	Value
Product Name	GSC Rabbit mAb
Remarks/Alias	SAMS; GSC
Species	Human
GeneID (Human)	145258
GeneID	145258
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human GSC (P56915)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P56915
Protein Weight	28kDa
Shipping	Ice bag