

## Anti-Gsdma Monoclonal Antibody For Wb, If-P, Elisa - Q96Qa5

Item Number: CM0006413

### Introduction

Rabbit monoclonal antibody against human GSDMA (Gasdermin-A). Validated for WB, IF-P, ELISA; cross-reacts with mouse. Suitable for pyroptosis research and epithelial immunity studies.

[Learn More](#)

| Parameter        | Value                     |
|------------------|---------------------------|
| Product Name     | GSDMA Rabbit mAb          |
| Remarks/Alias    | GSDM; FKSG9; GSDM1; GSDMA |
| Species          | Human                     |
| GeneID (human)   | 284110                    |
| GeneID           | 284110                    |
| Immunogen        | Recombinant protein       |
| Source           | Rabbit                    |
| Category         | Monoclonal Antibodies     |
| Application      | WB, IF-P, ELISA           |
| Cross Reactivity | Human, Mouse              |
| SWISS            | Q96QA5                    |
| Protein Weight   | 49kDa                     |
| Shipping         | Ice bag                   |