

Anti-Gsdmc Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Byg8

Item Number: CM0016185

Introduction

High-quality rabbit anti-GSDMC polyclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human and mouse GSDMC, a 58 kDa pore-forming pyroptosis protein. Suitable for studying gasdermin-C-mediated cell death pathways.

[Learn More](#)

Parameter	Value
Product Name	GSDMC Rabbit pAb
Remarks/Alias	MLZE; GSDMC
Species	Human
GeneID (Human)	56169
GeneID	56169
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human GSDMC (NP_113603.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BYG8
Protein Weight	58kDa
Shipping	Ice bag