


Anti-Gsdmd (Full Length+C Terminal) Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q9D8T2

Item Number: CM0004086



Introduction

Anti-GSDMD (Full Length+C terminal) rabbit monoclonal antibody validated for WB, IF/ICC, IHC-P, ELISA in mouse. Ideal for studying pyroptosis, inflammasome signaling, and gasdermin D pore formation.

[Learn More](#)

Parameter	Value
Product Name	GSDMD (Full Length+C terminal) Rabbit mAb
Remarks/Alias	DF5L; M2-4; Dfna5l; GsdmD-1; Gsdmdc1; 1810036L03Rik
Species	Mouse
Gene ID (Human)	79792
Gene ID	69146
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	Q9D8T2
Protein Weight	53 kDa
Shipping	Ice bag