


Anti-Gsdmd (Full Length+N Terminal) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P57764

Item Number: CM0022815



Introduction

Rabbit polyclonal antibody against human GSDMD (Full length+N terminal). Validated for WB, IF/ICC, ELISA with cross-reactivity to human, mouse, rat. Target is a pore-forming protein critical for pyroptosis and host defense.

[Learn More](#)

Parameter	Value
Product Name	GSDMD (Full length+N terminal) Rabbit pAb
Remarks/Alias	DF5L; DFNA5L; FKSG10; GSDMDC1; GSDMD (Full Length+N terminal)
Species	Human
GeneID (human)	79792
GeneID	79792
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P57764
Protein Weight	53kDa
Shipping	Ice bag