


Anti-Gsdmd (Full Length+N Terminal) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9D8T2

Item Number: CM0015118



Introduction

Rabbit polyclonal antibody targeting full-length and N-terminal GSDMD (Q9D8T2). Validated for WB, IF/ICC, ELISA in human, mouse, and rat. Ideal for pyroptosis and inflammasome research.

[Learn More](#)

Parameter	Value
Product Name	GSDMD (Full length+N terminal) Rabbit pAb
Remarks/Alias	DF5L; M2-4; Dfna5l; Gsdmd-1; Gsdmdc1; 1810036L03Rik; GSDMD (Full Length+N terminal)
Species	Mouse
Gene ID (Human)	79792
Gene ID	69146
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9D8T2
Protein Weight	53kDa
Shipping	Ice bag