


## Anti-Gsdmd Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P57764

Item Number: CM0024620



### Introduction

Rabbit polyclonal antibody targeting human GSDMD (Gasdermin-D), validated for WB, IHC-P, IF/ICC, IP, ELISA; cross-reacts with mouse, rat. For pyroptosis, inflammasome, and innate immunity research.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | GSDMD Rabbit pAb                                                                                                                        |
| Remarks/Alias    | DF5L; DFNA5L; FKSG10; GSDMDC1; GSDMD                                                                                                    |
| Species          | Human                                                                                                                                   |
| GeneID           | 79792                                                                                                                                   |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 346-484 of human GSDMD (NP_079012.3). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                                   |
| Application      | WB, IHC-P, IF/ICC, IP, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                       |
| SWISS            | P57764                                                                                                                                  |
| Protein Weight   | 53 kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                 |