

Anti-Gsk-3A/B Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P49840/P49841

Item Number: CM0006581



Introduction

Rabbit monoclonal antibody targeting human GSK-3 α and GSK-3 β . Validated in WB, IF/ICC, IF-P, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for diabetes, Alzheimer's, and Wnt signaling research.

[Learn More](#)

Parameter	Value
Product Name	GSK-3 α / β Rabbit mAb
Remarks/Alias	GSK3B; gsk-3 β ; GSK-3 α / β
Species	Human
GeneID (Human)	2931/2932
GeneID	2931/2932
Immunogen	Recombinant protein of human GSK-3 α / β
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49840/P49841
Protein Weight	46kDa/48kDa
Shipping	Ice bag