


Anti-Gsk-3A/B Monoclonal Antibody For Wb, Ip, Elisa - P49840 / P49841

Item Number: CM0004001



Introduction

Anti-GSK-3 α / β Rabbit Monoclonal Antibody validated for WB, IP, and ELISA in human, mouse, and rat. Ideal for studying Wnt signaling, glucose homeostasis, and Alzheimer disease pathways.

[Learn More](#)

Parameter	Value
Product Name	GSK-3 α / β Rabbit mAb
Remarks/Alias	GSK3B; gsk-3 β ; GSK-3 α / β
Species	Human
GeneID (Human)	2931/2932
GeneID	2931/2932
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49840/P49841
Protein Weight	46kDa/48kDa
Shipping	Ice bag