


Anti-Gsk-3A/B Rabbit Polyclonal Antibody For Ihc-P, Elisa - P49840 / P49841

Item Number: CM0021056



Introduction

Rabbit polyclonal antibody targeting human and mouse GSK-3 α / β (P49840/P49841) for IHC-P, ELISA applications. Useful for studies on Wnt signaling, glycogen metabolism, Alzheimer's disease, and autophagy.

[Learn More](#)

Parameter	Value
Product Name	GSK-3 α / β Rabbit pAb
Remarks/Alias	GSK3B; gsk-3 β ; GSK-3 α / β
Species	Human
Gene ID (Human)	2931/2932
Gene ID	2931/2932
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49840/P49841
Protein Weight	46kDa/48kDa
Shipping	Ice bag