


Anti-Gsk3A Polyclonal Antibody For Wb - P49840

Item Number: CM0018784



Introduction

Rabbit polyclonal antibody targeting human GSK3 α , validated for WB and ELISA. Cross-reacts with human, mouse, rat. Ideal for Wnt signaling, insulin resistance, and Alzheimer's disease research. ~51 kDa.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GSK3 α Rabbit pAb
Remarks/Alias	GSK3A; glycogen synthase kinase-3 alpha; 3 α
Species	Human
Gene ID (Human)	2931
Gene ID	2931
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 400-483 of human GSK3 α (NP_063937.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49840
Protein Weight	51kDa
Shipping	Ice bag