


Anti-Gsk3B Monoclonal Antibody For Flow Cytometry (Intracellular) - P49841

Item Number: CM0005589



Introduction

-conjugated rabbit monoclonal antibody against human GSK3 β for intracellular flow cytometry. Targets glycogen synthase kinase-3 beta; no cross-reactivity data provided. For research use in Wnt signaling, metabolism, and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human GSK3 β mAb
Remarks/Alias	GSK3B; glycogen synthase kinase-3 beta; GSK3 β
Species	Human
GeneID (human)	2932
GeneID	2932
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 321-420 of human GSK3 β (NP_001139628.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P49841
Protein Weight	47kDa
Shipping	Ice bag