


Anti-Gsk3B Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P49841

Item Number: CM0028930



Introduction

High-quality GSK3 β rabbit polyclonal antibody validated for WB, IF/ICC, IP, and ELISA applications. Reacts with human, mouse, and rat GSK3 β (47 kDa). Suitable for Wnt signaling, insulin resistance, and Alzheimer's disease research.

[Learn More](#)

Parameter	Value
Product Name	GSK3 β Rabbit pAb
Remarks/Alias	GSK3B; gsk-3 β ; GSK3 β
Species	Human
GeneID (Human)	2932
Immunogen	Synthetic peptide (aa 320-420 of human GSK3 β , NP_001139628.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49841
Protein Weight	47kDa
Shipping	Ice bag