


Anti-Gspt1 Rabbit Polyclonal Antibody For Wb, Elisa - P15170

Item Number: CM0019817



Introduction

Rabbit polyclonal antibody against human GSPT1 (ERF3A) for WB and ELISA. Detects human, mouse, and rat GSPT1. Relevant to translation termination and nonsense-mediated decay studies.

[Learn More](#)

Parameter	Value
Product Name	GSPT1 Rabbit pAb
Remarks/Alias	GST1; ETF3A; eRF3a; 551G9.2; GSPT1
Species	Human
Gene ID (Human)	2935
Gene ID	2935
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 138-420 of human GSPT1 (NP_001123478.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15170
Protein Weight	56kDa
Shipping	Ice bag