

Anti-Gsr Rabbit Monoclonal Antibody For Wb, Elisa - P00390

Item Number: CM0008205



Introduction

High-performance Anti-GSR rabbit monoclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat GSR. Ideal for studying glutathione metabolism, oxidative stress, and redox biology research.

[Learn More](#)

Parameter	Value
Product Name	GSR Rabbit mAb
Remarks/Alias	GR; GSRD; HEL-75; HEL-S-122m; GSR
Species	Human
Gene ID	2936
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 423-522 of human GSR (P00390).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00390
Protein Weight	56kDa
Shipping	Ice bag