


Anti-Gst3/Gstp1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P09211

Item Number: CM0003718



Introduction

Rabbit monoclonal antibody against human GSTP1 (GST3/GSTP1), validated for WB, IHC-P, ELISA, cross-reacts with mouse and rat. Detects ~23kDa band. Ideal for detoxification, ferroptosis, neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	GST3/GSTP1 Rabbit mAb
Remarks/Alias	PI; DFN7; GST3; GSTP; FAEES3; HEL-S-22; GST3/GSTP1
Species	Human
GeneID (Human)	2950
GeneID	2950
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09211
Protein Weight	23kDa
Shipping	Ice bag