


Anti-Gst3/Gstp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09211

Item Number: CM0024856



Introduction

Rabbit polyclonal antibody targeting human GST3/GSTP1, a key enzyme in glutathione conjugation and ferroptosis regulation. Suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reactive with mouse. Ideal for studies on oxidative stress and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	GST3/GSTP1 Rabbit pAb
Remarks/Alias	PI; DFN7; GST3; GSTP; FAEES3; HEL-S-22; GST3/GSTP1
Species	Human
Gene ID	2950
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human GST3/GST3/GSTP1 (NP_000843.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09211
Protein Weight	23kDa
Shipping	Ice bag