

Anti-Gsta1 Polyclonal Antibody For Wb, If/Icc, Elisa - P08263

Item Number: CM0018868



Introduction

High affinity polyclonal antibody against human GSTA1, validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for detoxification and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	GSTA1 Rabbit pAb
Remarks/Alias	GST2, GTH1, GSTA1-1, GST-epsilon, GSTA1
Species	Human
Gene ID (Human)	2938
Gene ID	2938
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-222 of human GSTA1 (NP_665683.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08263
Protein Weight	26kDa
Shipping	Ice bag