

Anti-Gsta1 Rabbit Polyclonal Antibody For Wb, Elisa - P08263

Item Number: CM0019564



Introduction

Rabbit polyclonal antibody targeting human GSTA1 (glutathione S-transferase A1). Validated for WB and ELISA, with cross-reactivity in rat. Ideal for detoxification and hormone biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	GSTA1 Rabbit pAb
Remarks/Alias	GST2; GTH1; GSTA1-1; GST-epsilon; GSTA1
Species	Human
GeneID (Human)	2938
GeneID	2938
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-222 of human GSTA1 (NP_665683.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P08263
Protein Weight	26kDa
Shipping	Ice bag