

Anti-Gstk1 Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Y2Q3

Item Number: CM0004419



Introduction

Recombinant rabbit monoclonal antibody against human glutathione S-transferase kappa 1 (GSTK1). Validated for Western blot, immunofluorescence, and ELISA with cross-reactivity to mouse. Suitable for research on peroxisomal transferase function.

[Learn More](#)

Parameter	Value
Product Name	GSTK1 Rabbit mAb
Remarks/Alias	GST; GST13; hGSTK1; GSTK1-1; GST13-13; GST 13-13; GSTK1
Species	Human
GeneID (Human)	373156
GeneID	373156
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human GSTK1 (Q9Y2Q3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y2Q3
Protein Weight	25kDa
Shipping	Ice bag