

Anti-Gstk1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y2Q3

Item Number: CM0021250



Introduction

Polyclonal antibody against human GSTK1, validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for studying glutathione metabolism and peroxisomal detoxification pathways.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GSTK1 Rabbit pAb
Remarks/Alias	GST; GST13; hGSTK1; GSTK1-1; GST13-13; GST 13-13; K1
Species	Human
Gene ID (Human)	373156
Gene ID	373156
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-226 of human GSTK1 (NP_057001.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2Q3
Protein Weight	25kDa
Shipping	Ice bag