

## Anti-Gstk1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y2Q3

Item Number: CM0026918



### Introduction

Anti-GSTK1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA in human, mouse, and rat samples. Detects the 25kDa peroxisomal glutathione S-transferase kappa 1, ideal for detoxification and subcellular localization studies.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | GSTK1 Rabbit pAb                                                                                                 |
| Remarks/Alias    | GST; GST13; hGSTK1; GSTK1-1; GST13-13; GST 13-13; GSTK1                                                          |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 373156                                                                                                           |
| GeneID           | 373156                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-226 of human GSTK1 (NP_057001.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q9Y2Q3                                                                                                           |
| Protein Weight   | 25kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |