

Anti-Gstm1 Polyclonal Antibody For Wb, Elisa - P09488

Item Number: CM0024190



Introduction

High-quality anti-GSTM1 polyclonal antibody validated for Western blot (WB) and ELISA. Detects human, mouse, and rat GSTM1. Ideal for studying glutathione metabolism and detoxification pathways.

[Learn More](#)

Parameter	Value
Product Name	GSTM1 Rabbit pAb
Remarks/Alias	MU; H-B; GST1; GTH4; GTM1; MU-1; GSTM1-1; GSTM1a-1a; GSTM1b-1b; GSTM1
Species	Human
Gene ID	2944
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-181 of human GSTM1 (NP_666533.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09488
Protein Weight	26kDa
Shipping	Ice bag