

Anti-Gstm2 Rabbit Polyclonal Antibody For Wb, Elisa - P28161

Item Number: CM0013985



Introduction

Rabbit polyclonal antibody targeting human GSTM2, a glutathione S-transferase involved in detoxification and lipid metabolism. Validated in WB and ELISA with cross-reactivity to mouse. Ideal for metabolic and toxicology research.

[Learn More](#)

Parameter	Value
Product Name	GSTM2 Rabbit pAb
Remarks/Alias	GST4; GSTM; GTHMUS; GSTM2-2; GSTM2
Species	Human
Gene ID (Human)	2946
Gene ID	2946
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of human GSTM2 (NP_000839.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P28161
Protein Weight	26kDa
Shipping	Ice bag