

Anti-Gstm3 Polyclonal Antibody For Wb, If/Icc, Elisa - P21266

Item Number: CM0028108



Introduction

Rabbit polyclonal antibody against human GSTM3 (GST class-mu 3). Suitable for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. For detoxification and xenobiotic metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	GSTM3 Rabbit pAb
Remarks/Alias	GST5; GSTB; GTM3; GSTM3-3; GSTM3TV2; hGSTM3-3; GSTM3
Species	Human
Gene ID	2947
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 101-225 of human GSTM3 (NP_000840.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21266
Protein Weight	27kDa
Shipping	Ice bag